

## OxiTop BOD Measurement Instrumentation

Pressure sensors for BOD testing and measurements.



- Simple operation
- Non-Toxic
- Stores 100 user-defined routine measurements
- User-friendly on screen guidance for easy operation

### Description

#### OxiTop IS 6, IS 12

Measurement using OxiTop® Biological Oxygen Demand (BOD) Instrumentation is based on pressure measurement in a closed system: microorganisms in the sample consume the oxygen and form CO<sub>2</sub>. This is absorbed by NaOH, creating a vacuum, which can be read directly as a measured value in mg/l BOD.

The sample volume being tested regulates the amount of oxygen available for a complete BOD measurement. BOD measurement ranges of up to 4,000 mg/l can be measured using different sample volumes.

The OxiTop® BOD instrumentation heads (green and yellow for differentiation of inflow/outflow) have an AutoTemp function: if the sample temperature is too cold, the start of BOD measurement is automatically delayed (by at least 1 hour) until a constant temperature has been reached.

Apart from the automatic storage of 5 measured values (1 value per day), further measured BOD values can be read at all times during or after the period of 5 days, which permits the tracking of check values or measurements over longer periods. The OxiTop BOD Instrumentation measurement system also offers simple operation, improved controllability, and non-toxicity.

### Specifications

*Measuring principle:* Respirometric (manometric)  
*Display:* 2-digit 7-segment LED, 10 mm  
*Measured parameter:* Digits  
*Measuring range:* 0 to 40 digits (+10 digits for overranging)  
*Accuracy:* ± 0.03 psi (± 2 hPa) or ± 1 digit  
*Admissible measuring temperature:* 68°F ± 1.8°F (20°C ± 1°C)  
*Admissible sample temperature (when being filled):* 59-68°F (15-20°C)  
*Power supply:* Lithium batteries (260 mAh), 2 x type CR 2430 (3V)  
*Power consumption:* Max. 25 mA (during measurement)  
*Battery run time:* ~ 3 years with normal use  
*Protection class:* 3 IEC 1010  
*Protection system:* IP 54 IEC 529  
**EMC**  
*Emissions:* EN 50081-1, FCC Class A  
*Immunity:* EN 50082-2, NAMUR recommendation  
*Climate class:* 2, VDI/VDE 3540  
*Ambient temperature:*  
*Storage:* -13°F to +149°F (-25°C to +65°C)  
*Operation:* 68°F ± 1.8°F (20°C ± 1°C)  
*Test mark:* CE  
*Dimensions:* H: 2.7in, dia. 2.8in (69 mm, dia. 70 mm)  
*Weight:* Approx. 3 oz. (85 g)

### Options and Accessories

#### 208214 OxiTop® IS 6

BOD Measurement Instrumentation  
 BOD measuring instrument for self-check with 6 measuring units, 115 VAC, 50/60 Hz power plug, and accessories.

#### 208215 OxiTop® IS 12

BOD Measurement Instrumentation  
 BOD measuring instrument for self-check with 12 measuring units, 115 VAC, 50/60 Hz power plug, and accessories.

#### 208433 OxiTop® Thermostat Box

Thermostat Box  
 OxiTop® BOD Instrumentation thermostat box, for up to 12 OxiTop® units and stirrer base, 115 VAC, 50 Hz.

Contact  
 Global Water  
 for all your  
 instrumentation  
 needs:

Water Level

Water Flow

Water Samplers

Water Quality

Weather

Remote Monitoring

Control